

Product Datasheet

RFC4 antibody (orb648689)

Description	Rabbit polyclonal antibody to RFC4
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	Recombinant full length Human RFC4.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~39 kDa
Uniprot ID	P35249
Dilution Range	WB1:500 - 1:2000, IHC-P: -P1:50 - 1:200, IF1:50 - 1:100IP1:50 - 1:200
Expiration Date	12 months from date of receipt.