

Product Datasheet

p40phox (phospho-Thr154) antibody (orb6491)

Description	p40phox (phospho-Thr154) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human p40phox around the phosphorylation site of Thr154 (PR(p-T)RK)
Target	Phospho-NCF4 (Thr154)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	Q15080
Expiration Date	12 months from date of receipt.