

Product Datasheet

SR-2C (P190) antibody (orb650649)

Description	Rabbit polyclonal antibody to SR-2C (P190)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic peptide, corresponding to amino acids 160-210 of Human SR-2C.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 55 kDa
Uniprot ID	P28335
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.