

Product Datasheet

ATP-Citrate Lyase antibody (orb651138)

Description	Rabbit monoclonal antibody to ATP-Citrate Lyase
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human ACLY(NP_001087.2).
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Form/Appearance	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	125KDa
Uniprot ID	P53396
Dilution Range	WB: 1:1000-1:2000 , , IF: 1:50-1:200, FC: 1:50-1:100
Expiration Date	12 months from date of receipt.