

Product Datasheet

HA-tag (26D11) antibody (orb651177)

Description	Mouse monoclonal antibody to HA-tag (26D11)
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant full length Human GPR87.
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Form/Appearance	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB: 1:5000-1:10000
Expiration Date	12 months from date of receipt.