

## Product Datasheet

### mCherry-Tag antibody (orb651333)

<b>Description</b>	Mouse monoclonal antibody to mCherry-Tag
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	All
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, IHC, WB
<b>Immunogen</b>	Synthetic peptide mCherry-Tag.
<b>Preservatives</b>	Mouse IgG1, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	WB: 1:5000-1:10000, IHC-P: 1:200, IF: 1:200
<b>Expiration Date</b>	12 months from date of receipt.