

Product Datasheet

RFP-Tag antibody (orb651334)

Description	Mouse monoclonal antibody to RFP-Tag
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide RFP-Tag.
Preservatives	Mouse IgG1, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB: 1:5000-1:10000, IHC-P: 1:200, IF: 1:200
Expiration Date	12 months from date of receipt.