

## Product Datasheet

### V5-tag antibody (orb651348)

<b>Description</b>	Mouse monoclonal antibody to V5-tag
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	All
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Synthetic peptide, corresponding to V5 tag.
<b>Preservatives</b>	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>MW</b>	~ 95 kDa
<b>Dilution Range</b>	WB: 1:5000-1:20000IP:1:100-1:1000
<b>Expiration Date</b>	12 months from date of receipt.