

Product Datasheet

MLST8 antibody (orb652072)

Description	Mouse monoclonal antibody to MLST8
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	MLST8
Form/Appearance	Mouse monoclonal antibody supplied in crude ascites with 0.09% (W/V) sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	35876
Uniprot ID	Q9BVC4
NCBI	NP_001186103.1 , NP_071767.3 , NP_001186104.1 , NP_001186102.1
Dilution Range	WB: 1:100-8000
Expiration Date	12 months from date of receipt.