

Product Datasheet

PLOD1 antibody (orb652403)

Description	Mouse monoclonal antibody to PLOD1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	PLOD1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG2b, κ
Clonality	Monoclonal
MW	83550
Uniprot ID	Q02809
Dilution Range	WB: 1:4000
Expiration Date	12 months from date of receipt.