

Product Datasheet

NECAB3 antibody (orb652553)

Description	Rabbit polyclonal antibody to NECAB3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q96P71
Dilution Range	IF: 1:100
Expiration Date	12 months from date of receipt.