



Product Datasheet

MMP10 antibody (orb652596)


Description	Rabbit polyclonal antibody to MMP10
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P09238
Dilution Range	IF: 1:100, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cel...



Western blot analysis of A549 cell line,...



Immunohistochemistry analysis of paraffi...