

## Product Datasheet

### KDEL2 antibody (orb652614)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Rabbit polyclonal antibody to KDEL2  |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human, Mouse, Rat  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | IF, WB   |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| <b>Note</b>                | For research use only  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <a href="#">P33947</a>   |
| <b>Dilution Range</b>      | IF: 1:100  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |