


## Product Datasheet

### Cytochrome P450 27A1 Antibody (orb652684)

|                            |                                                                                                                                        |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to CP27A                                                                                                    |
| <b>Species/Host</b>        | Rabbit                                                                                                                                 |
| <b>Reactivity</b>          | Human                                                                                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                                                                                           |
| <b>Tested Applications</b> | WB                                                                                                                                     |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| <b>Note</b>                | For research use only                                                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                                                                             |
| <b>Uniprot ID</b>          | <b>Q02318</b>                                                                                                                          |
| <b>Dilution Range</b>      | WB: 1:500                                                                                                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                        |



Western blot analysis of lysate from human liver tissue lysate, using Cytochrome P450 27A1 Antibody. Diluted at 1:1000. A goat anti-rabbit IgG H&L (HRP) at 1:5000 dilution was used as the secondary antibody. Lysate at 35 ug.

#### **Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### **Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)