

Product Datasheet

TAB1 antibody (orb653416)

Description	Mouse monoclonal antibody to TAB1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q15750
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.