

Product Datasheet

CKMT1A antibody (orb653429)

Description	Mouse monoclonal antibody to CKMT1A
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P12532
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.