


Product Datasheet

Beclin 1 Antibody (orb653437)

Description	Mouse monoclonal antibody to Beclin 1
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q14457
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot detection of Beclin 1 in HeLa, 3T3, C6, COS7, CHO-K1, Rat brain and Jurkat cell lysates using Beclin 1 mouse mAb (1:1000 diluted). Predicted band size: 52KDa. Observed band size: 55KDa.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)