

Product Datasheet

EGFR antibody (orb653503)

Description	Mouse monoclonal antibody to EGFR
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P00533
Dilution Range	IP: 1:500
Expiration Date	12 months from date of receipt.