



Product Datasheet

S100a10 antibody (orb653589)


Description	Mouse monoclonal antibody to S100a10
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P08207
Dilution Range	ICC: 1:200, IP: 1:100
Expiration Date	12 months from date of receipt.



Western blot detection of S100A10 in C6,...



Immunocytochemistry staining of HeLa cel...



Immunoprecipitation analysis of HeLa cel...