


Product Datasheet

COVID-19 Nucleocapsid protein antibody (orb653664)

Description	Recombinant human monoclonal antibody recognizes COVID-19 Nucleocapsid protein.
Species/Host	Human
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Target	N
Preservatives	1 mg/ml. IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.4.
Form/Appearance	1 mg/ml. IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Human Cells
Uniprot ID	P0DTC9
Entrez	43740575
Dilution Range	ELISA: Use at an assay dependent dilution.
Expiration Date	12 months from date of receipt.



Indirect
ELISA
analysis of
COVID-19
Nucl...