


Product Datasheet

COVID-19 Nucleocapsid Protein scFv Antibody (orb653665)

Description	Recombinant monoclonal antibody is produced by E.coli expression system with a 6His, Flag tag at the...
Species/Host	E. coli
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Target	N
Preservatives	1 mg/ml. IgG1. Liquid in 0.02M Phosphate Buffered Saline, pH 7.4.
Form/Appearance	1 mg/ml. IgG1. Liquid in 0.02M Phosphate Buffered Saline, pH 7.4.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	E.coli
Uniprot ID	PODTC9
Entrez	43740575
Dilution Range	ELISA: Use at an assay dependent dilution.
Expiration Date	12 months from date of receipt.



Indirect
ELISA
analysis of
COVID-19
Nucl...