

Product Datasheet

Spike glycoprotein antibody (orb653753)

Description	Mouse monoclonal antibody to Spike glycoprotein
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, GICA
Immunogen	Recombinant Human Novel Coronavirus Spike glycoprotein (S) (16-685aa)
Target	S
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Monoclonal mouse (variable region) / human (kappa / IgG1 constant) chimeric antibody
Clonality	Monoclonal
Uniprot ID	P0DTC2
Dilution Range	ELISA:1:10000-1:50000:1:500-1:25000
Expiration Date	12 months from date of receipt.