

Product Datasheet

SARS-CoV-2 Spike RBD N antibody (orb653755)

Description	Mouse monoclonal antibody to SARS-CoV-2 Spike RBD N
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, GICA, NeA
Immunogen	Recombinant Human Novel Coronavirus Spike glycoprotein(S) (319-541aa)
Target	S
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	VHH fusion with human IgG1 Fc
Clonality	Monoclonal
Uniprot ID	P0DTC2
Dilution Range	ELISA:1:10000-1:100000:1:10000-1:40000
Expiration Date	12 months from date of receipt.