

Product Datasheet

Nucleoprotein antibody (orb653760)

Description	Mouse monoclonal antibody to Nucleoprotein
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, GICA, WB
Immunogen	Recombinant Human Novel Coronavirus Nucleoprotein (N) (1-419aa)
Target	N
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse scFv fusion with human IgG1 Fc
Clonality	Monoclonal
Uniprot ID	P0DTC9
Dilution Range	ELISA:1:10000-1:50000, WB:1:500-1:5000:1:1000-1:10000
Expiration Date	12 months from date of receipt.