

Product Datasheet

Nucleoprotein antibody (HRP) (orb653761)

Description	Mouse monoclonal antibody against Nucleoprotein conjugated to HRP
Reactivity	Virus
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Novel Coronavirus Nucleoprotein (N) (1-419aa)
Target	N
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Mouse scFv fusion with human IgG1 Fc
Clonality	Monoclonal
Uniprot ID	P0DTC9
Dilution Range	Recommended dilution: ELISA:1:10000-1:50000
Expiration Date	12 months from date of receipt.