

Product Datasheet

HighFidelity ATTO594 PCR Labeling Kit (orb653769)

Description HighFidelity ATTO594 PCR Labeling Kit is designed to produce randomly ATTO594-modified DNA probes b...

Storage store at -20 °C

Note For research use only

Application notes **Spectroscopic Propertie:** λ_{exc} 602 nm, λ_{em} 626 nm, ϵ 120.0 L mmol⁻¹ cm⁻¹ (Tris-HCl pH 7.5).

SMILES C(/C=C/c1cn(c(=O)[nH]c1=O)[C@H]1C[C@@H]([C@H](O1)COP(=O)(O)OP(=O)(O)OP(=O)(O)O)NC(=O)CCCCNC(=O)CCCCNC(=O)CCCN(C(=O)c1c(ccc1)C1=c2c(Oc3c1cc1c(c3)N(C(C=C1CS(=O)(=O)O)(C)C)CC)cc1=[N+](C(C=C(c1c2)CS(=O)(=O)[O-])(C)C)CC)C

Expiration Date 12 months from date of receipt.