

## Product Datasheet

# HighFidelity ATTO643 PCR Labeling Kit (orb653770)

<b>Description</b>	HighFidelity ATTO643 PCR Labeling Kit is designed to produce randomly ATTO643-modified DNA probes b...
<b>Storage</b>	store at -20 °C
<b>Note</b>	For research use only
<b>Application notes</b>	<b>Spectroscopic Propertie:</b> $\lambda_{exc}$ 643 nm, $\lambda_{em}$ 665 nm, $\epsilon$ 150.0 L mmol <sup>-1</sup> cm <sup>-1</sup> (Tris-HCl pH 7.5).
<b>SMILES</b>	<chem>C(/C=C/c1cn(c(=O)[nH]c1=O)[C@H]1C[C@@H]([C@H](O1)COP(=O)(O)OP(=O)(O)OP(=O)(O)O)O)N*</chem>
<b>Expiration Date</b>	12 months from date of receipt.