

Product Datasheet

Goat IgG Protein (orb653870)

Description	goat immunoglobulin G
Tested Applications	SDS-PAGE
Preservatives	Phosphate buffer (0.1 M, pH 7.4) and 0.001% sodium azide
Concentration	10 mg/ml
Storage	Store at -80°C for long-term storage.
Note	For research use only
Application notes	The protein solution should be gently mixed before use. Avoid freeze/thaw cycles.
Purity	> 95% by SDS-PAGE
Expiration Date	6 months from date of receipt.