

Product Datasheet

His Tag antibody (orb653872)

Description	Mouse monoclonal antibody to Polyhistidine Tag
Conjugation	Unconjugated
Preservatives	1mg/ml in PBS (pH 7.4) 0.02% Sodium Azide & 10% Glycerol
Form/Appearance	Sterile Filtered solution.
Note	For research use only
Application notes	Applications: ELISA and W.B. The researcher should titrate to obtain optimal results
Isotype	Kappa IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.