

## Product Datasheet

### NMDAR2B (phospho-Tyr1336) antibody (orb6541)

<b>Description</b>	NMDAR2B (phospho-Tyr1336) antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	KLH conjugated synthesised phosphopeptide derived from human NMDAR2B around the phosphorylation site of Tyr1336 (SP(p-Y)AH)
<b>Target</b>	Phospho-NMDAR2B (Tyr1336)
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	164 kDa
<b>CAS Number</b>	Glutamate receptor ionotropic, NMDA 2B
<b>Uniprot ID</b>	<a href="#">Q13224</a>
<b>Dilution Range</b>	WB=1:500-2000
<b>Expiration Date</b>	12 months from date of receipt.