
Product Datasheet

Hexokinase 1/HK1 Antibody (monoclonal, 2I4) (orb654273)

Description

Hexokinase 1/HK1 Antibody (monoclonal, 2I4)

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

FC, ICC, IF, IHC, IHC-Fr, WB

Immunogen

E.coli-derived human Hexokinase 1/HK1 recombinant protein (Position: D17-R323).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Mouse IgG2b

Clonality


Monoclonal

MW


120 kDa

Uniprot ID
[P19367](#)
Expiration Date


12 months from date of receipt.



Western blot analysis of Hexokinase 1/HK...



IHC analysis of Hexokinase 1/HK1 using a...



IHC analysis of Hexokinase 1/HK1 using a...