
Product Datasheet

E3 SUMO-protein ligase PIAS1 PIAS1 Antibody (orb654360)

Description

E3 SUMO-protein ligase PIAS1 PIAS1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

ELISA, FC, ICC, IF, WB

Immunogen

E.coli-derived human PIAS1 recombinant protein (Position: Q96-K300).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality


Polyclonal

MW


72 kDa

Uniprot ID
[O75925](#)
Expiration Date


12 months from date of receipt.

Western blot analysis of PIAS1 using anti-PIAS1 antibody.



IF analysis of PIAS1 using anti-PIAS1 antibody.



Flow Cytometry analysis of PC-3 cells using anti-PIAS1 antibody.