

Product Datasheet

Biotin Conjugated AffiniPure Donkey Anti-Mouse IgM μ Chain (orb654483)

Description	0
Species/Host	Donkey
Conjugation	Biotin
Tested Applications	ELISA, IHC
Immunogen	Mouse IgM
Form/Appearance	Concentrated, Liquid
Concentration	1mg/ml
Storage	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Immunohistochemistry: 5-10 μ g/mL ELISA 0.05-0.2 μ g/mL. This antibody must be used with various kinds of enzymes or fluorescence conjugated avidin (streptavidin) in IHC and ELISA. Please centrifuge the vial before use. Dilute with neutral PBS or TBS. Prepare working solution on day of use.
Isotype	IgM
Expiration Date	12 months from date of receipt.