

## Product Datasheet

# HRP Conjugated AffiniPure Donkey Anti-Rabbit IgG(H+L) (orb654484)

<b>Description</b>	This peroxidase conjugate (HRP) conjugated antibody is specific for Rabbit IgG and shows no cross-reactivity with Bovine, Chicken, Goat, Guinea Pig, Syrian Hamster, Horse, Human, Mouse, Rat, and Sheep Serum Proteins.
<b>Species/Host</b>	Donkey
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Rabbit IgG (whole molecular)
<b>Form/Appearance</b>	Concentrated, Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
<b>Note</b>	For research use only
<b>Application notes</b>	Western Blotting: 0.005-1µg/ml (ECL detection) ELISA: 0.01-0.2µg/ml (TMB detection)
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)