

## Product Datasheet

### BMS-1001 (orb654513)

<b>Description</b>	A potent PD-1/PD-L1 interaction inhibitor with IC50 of 2.25 nM in a homogenous time-resolved fluores...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{35}H_{35}ClN_2O_7$
<b>MW</b>	631.122
<b>CAS Number</b>	[2113650-03-4]
<b>SMILES</b>	<chem>CC1=CC(=C(C=C1OCC2=C(C(=CC=C2)C3=CC4=C(C=C3)OCCO4)C)OC5=CC(=CC=C5)C#N)CN[C@H](CO)C(=O)O</chem>
<b>Expiration Date</b>	12 months from date of receipt.