

Product Datasheet

AUT-1 (orb654534)

Description	AUT-1 is a novel specific, cell permeant modulator of Kv3 channels with EC50 of 4.7 and 4.9 uM for K...
Note	For research use only
Formula	$C_{18}H_{19}N_3O_4$
MW	341.3611643
CAS Number	[1311136-84-1]
SMILES	<chem>C1(N(C2C=NC(OC3C=C(OC)C(C)=CC=3)=CC=2)C(=O)[C@H](CC)N1)=O</chem>
Expiration Date	12 months from date of receipt.