

Product Datasheet

Ustekinumab biosimilar antibody (orb654536)

Description	Ustekinumab is human monoclonal antibody to a polypeptide found on interleukin-12 and -23 that is us...
Conjugation	Unconjugated
Target	Ustekinumab
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
CAS Number	[815610-63-0]
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.