

## Product Datasheet

### Trastuzumab biosimilar antibody (orb654537)

<b>Description</b>	Trastuzumab (anti-human HER2) is a humanized, recombinant monoclonal antibody that binds to the extr...
<b>Endotoxins</b>	□1 EU/ml
<b>Conjugation</b>	Unconjugated
<b>Target</b>	Trastuzumab
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Bacterial Endotoxins Test: □1 EU/ml Exogenous Residual DNA: □1 pg/mg Residual protein A:□1 ng/mg
<b>Clonality</b>	Recombinant
<b>Purity</b>	Non-reduced CE-SDS: 97.5% SEC-HPLC: 99.7%
<b>CAS Number</b>	180288-69-1
<b>Hazard Information</b>	Non-Toxic
<b>Dilution Range</b>	Compared with standard, the range of biological activity is 110%
<b>Expiration Date</b>	12 months from date of receipt.