

Product Datasheet

Tocilizumab biosimilar antibody (orb654538)

| | |
|---------------------------|---|
| Description | Tocilizumab (anti-IL-6R) is a humanized monoclonal antibody that binds to the interleukin-6 receptor... |
| Conjugation | Unconjugated |
| Target | Tocilizumab |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Monoclonal |
| CAS Number | [375823-41-9] |
| Hazard Information | Non-Toxic |
| Expiration Date | 12 months from date of receipt. |