

Product Datasheet

Secukinumab biosimilar antibody (orb654539)

Description	Secukinumab is a fully human immunoglobulin G1-kappa monoclonal antibody that directly inhibits inte...
Conjugation	Unconjugated
Target	Secukinumab
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
CAS Number	[875356-43-7]
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.