

Product Datasheet

Pembrolizumab biosimilar antibody (orb654542)

Description	Pembrolizumab (anti-PD-1) is a potent, highly selective, fully humanized immunoglobulin (Ig) G4-kapp...
Conjugation	Unconjugated
Target	Pembrolizumab
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
CAS Number	[1374853-91-4]
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.