

Product Datasheet

Omalizumab biosimilar antibody (orb654543)

Description	Omalizumab is a monoclonal antibody targeting the high-affinity receptor binding site on human immun...
Conjugation	Unconjugated
Target	Omalizumab
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
CAS Number	[242138-07-4]
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.