

Product Datasheet

Nivolumab biosimilar antibody (orb654544)

Description Nivolumab anti-PD-1) is a genetically engineered fully human immunoglobulin (Ig) G4 monoclonal antibo...

Conjugation Unconjugated

Target Nivolumab

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Clonality Monoclonal

CAS Number [946414-94-4]

Hazard Information Non-Toxic

Expiration Date 12 months from date of receipt.