

Product Datasheet

Ipilimumab biosimilar antibody (orb654545)

Description	Ipilimumab (anti-CTLA-4) is an immunomodulatory monoclonal antibody directed against the cell surfac...
Conjugation	Unconjugated
Target	Ipilimumab
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
CAS Number	[477202-00-9]
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.