

Product Datasheet

Eculizumab biosimilar antibody (orb654547)

Description	Eculizumab is a long-acting humanized monoclonal antibody targeted against complement C5 and it is u...
Conjugation	Unconjugated
Target	Eculizumab
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
CAS Number	[219685-50-4]
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.