

Product Datasheet

Bevacizumab biosimilar antibody (orb654552)

Description	Bevacizumab (anti-VEGF) is a humanized anti-VEGF monoclonal antibody which binds to and neutralizes ...
Conjugation	Unconjugated
Target	Bevacizumab
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
CAS Number	[216974-75-3]
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.