

Product Datasheet

Alirocumab biosimilar antibody (orb654554)

Description	Alirocumab is a human monoclonal antibody that binds to proprotein convertase subtilisin kexin type ...
Conjugation	Unconjugated
Target	Alirocumab
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
CAS Number	[1245916-14-6]
SMILES	OC(CO)COP(O)(O)=O.[Na].[Na]
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.