

## Product Datasheet

### Aflibercept biosimilar antibody (orb654555)

<b>Description</b>	Aflibercept (VEGF-Trap) is a recombinant chimeric soluble receptor containing structural elements fr...
<b>Conjugation</b>	Unconjugated
<b>Target</b>	Aflibercept
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>CAS Number</b>	[862111-32-8]
<b>Hazard Information</b>	Non-Toxic
<b>Expiration Date</b>	12 months from date of receipt.