

Product Datasheet

Adalimumab biosimilar antibody (orb654556)

Description	Adalimumab (anti-TNF-alpha) is the first fully human, recombinant IgG1 monoclonal antibody that spec...
Conjugation	Unconjugated
Target	Adalimumab
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Recombinant
Purity	SEC-HPLC: 99.5%, Non-reduced CE-SDS: 98.7%
CAS Number	331731-18-1
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.