

Product Datasheet

BI4020 (orb654577)

Description	BI-4020 is a highly active, noncovalent EGFR inhibitor against EGFR mutations T790M and C797S with I...
Note	For research use only
Formula	$C_{30}H_{38}N_8O_2$
MW	542.675
SMILES	<chem>O=C1N=C2N(C3C=C(CN4CCN(C)CC4)C=CC=3N2)C[C@H](C)CCCOC2N(N=CC=2C2=CC1=CC(=N2)C)C t:2 </chem>
Expiration Date	12 months from date of receipt.